

# Ballancer® PLATINUM



Ballancer® lymphatic massage:  
The perfect addition to your body contouring treatments!





Superior aesthetic and  
wellbeing experience



Versatile, innovative  
programs



Simple to use:  
no training or specialized  
therapist needed



Intuitive touch screen interface to easily create unique massage programs

## Unique Features:

- Operates completely hands-free!
- Touch screen for easy programming
- Choose from preprogrammed treatment modes or create your own program
- Save and manage client lists and individual treatments
- Option of cancelling any chamber- full control of shorter cycles allows including more treatments



Ballancer® Pro Jacket



Ballancer® Pro Pants

## Treatment Modes:

Ballancer® Platinum preprogrammed treatment modes include:

- Cellulite programs for legs
- Intense program for arms
- Aesthetic regeneration
- Circulatory enhancement and muscular relaxation

Or, create your own massage sequences by combining Ballancer® lymphatic drainage modes:

- **Pretreatment** – primes the torso to maximize treatment results
- **Wave** – long, deep, and slow peristaltic massage
- **Ballancer®** – rapid, superficial, and fast caresses
- **Intense** – deep, directional lymphatic-drainage massage that starts at the foot or hand and proceeds to the torso, and then releases.

# Ballancer® Product Line

## How Ballancer® Works

Ballancer® uses special inflating compression garments to apply gentle or firm massage. Each Ballancer® garment contains 24 individual overlapping air chambers to apply fluent compression strokes that target each centimeter of the treated area. This massage motion, with its special direction of flow, increases circulation and gently but thoroughly decongests the problematic areas of hips, thighs, buttocks, and upper arms.

**The Ballancer® action accelerates the removal of waste products from fat metabolism, reducing the appearance of cellulite (“orange peel skin”), improving skin tone, and reducing volume in areas affected by excess fluid retention. Its gentle, rhythmic action promotes relaxation.**

## Unique Patented Garments

Our innovative, unique, and patent-protected garments are famous worldwide for their high quality, optimal treatment, and comfort.

- Treat entire upper or lower body
- Patent-protected fully overlapping chambers: unique system for smooth and effective treatment without gaps.
- More compression chambers: 24 compression chambers for efficient, fluent, comfortable massage
- Ballancer® Pro Pants can accommodate most sizes due to special design!
- Sturdy and durable fabrics that are easy to maintain and keep clean
- Advanced features for ultimate comfort and flexibility

## How Ballancer® Affects the Body



### CIRCULATORY SYSTEM

The Ballancer® peristaltic massage enhances blood circulation in the skin and deep layers under it. This brings fresher oxygen to skin cells and assists in speeding elimination of toxins.



### DIGESTIVE SYSTEM

Ballancer® compresses and massages the abdomen, relaxing abdominal muscles. The massage action is a natural bowel stimulant and can be especially helpful for clients who suffer from slow digestive system.



### LYMPHATIC SYSTEM

Ballancer® increases the lymphatic flow. The lymphatic system is responsible for removal of waste products from the areas under the skin, and Ballancer® massage speeds their elimination and fat discharge.



### MUSCLES

The Ballancer® peristaltic massage helps remove toxins from exercise and physical activity, as well as soothing, relaxing, and decongesting the muscles, easing aches and pains.



### SKIN

The Ballancer® massage increases the flow of freshly-oxygenated blood to skin layers. Skin will look firm and toned and will feel smoother and softer.

## Ordering Information

Part No.	Model
L10004CP	Ballancer® Platinum 1222
L4MP4136C	Ballancer® Pro Pants - Hip Hose Exit
L4MP4130C	Ballancer® Pro Pants - Heel Hose Exit
L4MP0001	Ballancer® Pro Disposable Pants
L4Z5BP03	Ballancer® Pro Pants - Expanders
L4MJ4120	Ballancer® Pro Jacket
L4Z54121	Ballancer® Pro Jacket - Expanders set Contains three expanders (1 back expander, 2 arm expanders)
L4MJ0001	Ballancer® Pro Disposable White Kimono